

T-JOE BRACE™ CARE AND MAINTENANCE



The T-Joe Brace is designed to allow freedom of movement and comfort. The brace uses lightweight materials to allow for the active lifestyle of kids and young adults involved in sports, gym class or just plain having fun!

The Chest Compression Plate™ of your brace is padded with orthopedic felt. The Back Support Pad™ uses our own design of padding called a bi-laminate and covered with our new Back Pad Covers. It is recommended that both be cleaned weekly. In your bathroom sink, mix warm water and antibacterial hand soap. Whip up the water with your fingers to form suds. Using a soft fingernail brush, scoop up some of the suds and gently scrub each pad. Repeat this step a few times. Next using a terry cloth towel pat them dry. Allow a few extra minutes to air dry before putting the brace back on. A blow dryer may also be used.

The Side Bar Sleeves that go over the bars on the left and right sides may be hand washed or put in a washing machine. **Do not over-dry as they may shrink.** Additional back pad covers and side bar sleeves are available from the T-Joe Pectus Store at www.PectusServices.com

The T-Joe Pectus Brace has been designed for maximum comfort which requires weekly maintenance. Each Monday gently hand tighten or “snug-up” the screws on your brace. To help you remember this we have termed the phrase *Maintenance Monday*. This means every Monday you will get a screwdriver and tighten up all the screws on your T-Joe Brace.

If any problems should arise contact our office for assistance.

THANK YOU FOR CHOOSING THE T-JOE PECTUS BRACING SYSTEM!



PECTUS SERVICES

PROVIDING PECTUS CARE SINCE 1996

877-PECTUS-6 (877-732-8876)

www.PectusServices.com