

CARE AND MAINTENANCE OF THE T-JOE PECTUS BRACING SYSTEM



The T-Joe brace is designed to allow freedom of movement and comfort. The brace uses lightweight materials to allow for an active lifestyle of fun playing sports, gym class at school. The padding of the brace can get dirty from sweat and skin oils we recommend they be cleaned as follows.

CHEST AND BACK PADS

For the Chest Pad and the Back Pad of your T-Joe Pectus Bracing System we recommend that they be cleaned weekly. To clean these pads, first set the stopper in your sink. Put in some warm water and a few drops of antibacterial hand soap. Then using your fingers 'slosh' the water to form suds. Use a small soft brush scoop up some of the suds and gently scrub the beige colored surface of the pads. Repeat this step a few times using clear water on the brush for the last one. Use an absorbent terry cloth towel to press and pat the pad dry. A blow dryer may be used on a cool setting for additional drying.

COVERS

The covers of your T-Joe Brace can simply be tossed into the washing machine. It is recommended that a fresh, clean set of covers be put on every 3 days or sooner if sweaty from sports or hot weather. When drying the covers use a cooler setting and shorter drying time to avoid shrinking. Additional sets of covers can be purchased through the T-Joe Store on our website: PectusServices.com.

If you should have any further questions or a problem should occur such as mechanical, hardware or padding contact your brace provider.

